

# Meters for X-Ray Service and QA/QC

## **QUART didoNEO**



## QUART *didoNEO* x-ray meter series

QUART GmbH  
Munich/Zorneding

Introducing the <i>QUART didoNEO</i>	3
<i>QUART didoNEO</i> Special Features	4
QUART Smart Accessories	5
History of QUART innovations	6

**dido** [/'dɪɪ.dɔ:/] (abbr.) shortened form of “diagnostic dosimeter”  
traditional ID for QUART meters

**NEO** [/'neɪ.ɔ:/] (ancient Greek) new, recent, current

**The QUART didoNEO  
line of meters introduce a new  
approach to the market of diagnostic x-ray measurement**

- The system features the **smallest** and **thinnest** multi-parameter detector available.
- It has the **lightest** multi-functional base unit ever designed in our industry.
- The user can access a **waveform preview** on the unit's display in the field – without the need to access a PC or laptop.
- Up to **10.000 exposures** can be **stored** in full for future reference or reporting.
- Together with the system, **several system options** are available providing **functional enhancement**.
- Technical **upgrades** will future proof the meters to cater for changes in user requirements.



Actual Size

## QUART didoNEO Special Features

### Compact Design Concept



Base unit and detector optimised in size and weight.

Most compact x-ray meter sensor in the industry.

Measurement behind scatter radiation grid for equipment attenuation factor / very low influence on fluoroscopy AEC.

Efficient vertical positioning (dental OPG)\*.

### Instant HVL and Total Filtration Measurement



Integrated direct-HVL measurement.

Instant total filtration measurement.

Both features are integrated and not charged extra.

### Exposure Data Storage Capacity



- Integrated long-time memory storage capacity.
- Storage of up to 10.000 measurements (depending on exposure event).

### Wave Form Preview



Tilt of the device provides exposure waveform.

Zoom-in and out for detailed visual analysis of the exposure graph.

### QUART NEOtec - Data Management and Reporting Software Tool



Transfer all measured data from didoNEO into the NEOtec software module for data management and analysis. Data base enables customer management, x-ray equipment data organisation and long-term data storage. Further Microsoft Excel export enables customised reporting and hardcopy printing.

### Accurate Dental OPG Measurement (The Dose-Width Product)



Wide sensor-surface design for accurate measurement at dental panoramic x-ray equipment. didoNEO automatically measures the dose-width product (DWP) of the OPG beam\*\*. Method accuracy as featured in didoNEO and didoEASY detectors officially verified by the University of Glasgow\*.

### Full Automatic Compensation



The didoNEO system is equipped with a multi-functional solid state detector. Measures almost any beam quality used in

diagnostic x-ray technology without preset procedure.

Minimal user interaction required.

### Exchangeable Detectors



To future-proof the meter for new requirements of our users additional detectors will be available for the didoNEO system (Release pending). These will cover:

- Mammography
- mA/mAs
- DLP in CT
- Survey Dose Measurements.



\* S A Mitchell and C J Martin, Comparison of ionisation chamber and semiconductor detector devices for measurement of the dose-width product for panoramic dental units, J. Radiol. Prot. 33 321 (2013)

\*\* In 1999 the National Radiological Protection Board (NRPB) recommended the introduction of dose-width product (DWP) for the measurement of patient dose in panoramic dental radiology. The DWP has further been recommended as a dose reference: Napier ID. Reference doses for dental radiography. Br Dent J (1999) 186: 392–6.

## QUART Smart Accessories

### QUART Bridge Holder

- ✓ Developed for easy and quick vertical positioning of our detectors.
- Holder provides superb fix on any smooth surface.
- Detector's padded anti-slip back prevents skidding on surfaces.
- Optional accessory.

### Base Unit Support

- ✓ Base stand for display readability from a distance.
- Enables tracking of displayed dose rate from a distance while measurement is running (e.g. in fluoroscopy applications).
- Quick set-up.
- Base stand is included in delivery.

### Hardshell Transport Case

- ✓ Robustly constructed case for all-time protection of your didoNEO - even when conditions are challenging.
- Bespoke foam insert designed to accept all system components plus accessories.
- Case is included in delivery.

### Hang Case

- ✓ Hang case provides safe storage for the meter where placement options are limited
- Viewing window provides reading the data while measurement is running.
- Optional accessory.



Hang Case



Bridge Holder



Robust Hardshell Transport Case



Customized Design Base Unit Stand

## History of QUART innovations

### QUART Integrator I

For diagnostic x-ray measurements, the Integrator I was the major step forward in the x-ray meter industry. The **world-wide first application of solid-state technology** in QUART's x-ray detectors, instead of utilising ionisation chambers, changed the characteristics of x-ray test equipment considerably. Once bulky and heavy, meters could now be easily handled and transported.



1975



1983

### QUART dido

The Integrator II, one year later (1984) renamed to QUART dido, became the conceptual base of all future QUART meters. The dido was the **first PTB\* approved diagnostic dosimeter of its kind** - featuring solid-state detector technology. Equipped with a slanted display for better reading from a distance or from above, the QUART dido set a new standard for shape and design of x-ray QA testing devices already at that time.

### QUART SPvario

The introduction of x-ray image quality assurance standards in the 1980s required the development of comprehensive but easy-to-use tools. Thanks to QUART's expertise and understanding, a **new generation of technical x-ray QA phantoms** was available from the beginning of this process.



1985



1988

### QUART dido/time | QUART R6Vi

Some time after its launch, the QUART R6Vi/time was further developed to become the **first sandwich/ double dosimeter**. It could be used to measure dose before and after patient equivalent filtration in constancy tests – in only one x-ray exposure.

### The DAVID System

The DAVID system for the first time featured a compact laptop computer as a **waveform analysis tool** to replace oscilloscopes previously used. Designed as a sophisticated and complete measurement system, it was the perfect tool for service experts and state radiation inspectors. After its launch it gained a reputation for causing a "toolbox revolution" in x-ray quality control. The system name DAVID transcribes as „**Digital Analyser for High-Voltage, Inherent Filtration and Dose Rate**“.



1992

\* The German PTB (Physikalisch-Technische Bundesanstalt) is the National Admission Authority to issue Type Approvals for measuring equipment. In Europe, the type approval usually is obligatory for meters used in commissioning tests.

<sup>1</sup> Deutsches Institut für Normung

<sup>2</sup> International Electrotechnical Commission

### Dental Test Phantom

The development of image quality test phantoms was launched quickly after the dental x-ray industry started inquiring for respective solutions. Until today, the general design features of QUART's dentFS and dent/digitest phantoms are still part of the **DIN<sup>1</sup> & IEC<sup>2</sup> standard for dental image quality control.**

1992



1993

### QUART SPdigi

The transformation from conventional to digital x-ray required a new approach on technology for x-ray QA. The **QUART-developed and standardised SPdigi phantom** has been incorporating test objects to specifically assess image quality of digital x-ray equipment.

### Digital Subtraction Angiography Phantom

The introduction of QUART's DSA phantom featuring longitudinal sliding technology has enabled a **precise** way to assess the imaging quality of subtraction angiography equipment. The method is still up-to-date and widely used. With a dedicated protocol, the test images can be evaluated using **modern digital methods.**

1996



2004 – 2005

### QUART dido2000K and dido2100K

The *dido2000K/2100K* series dosimeters are **all-in-one instruments** that incorporate kV and dose rate measurement. With their optional feature to output data via an USB interface, they further enabled waveform analysis and protocol print-outs.

### QUART ConeBeam CT Phantom and Software

The combination of phantom and software to analyse purely digital parameters, such as MTF and CNR, introduced a **whole new concept** into x-ray QA/QC. The software automatically evaluates phantom images and thus objectively assesses the imaging performance of CBCT x-ray systems.

2008



2012

### QUART didoNEO Series

The new didoNEO continues to advance the role that **genuine technology plays in latest measurement applications.** It expands user capabilities, maximises efficiency, increases flexibility and thus improves quality control and service workflow. The system was first presented in 2012.



**We help to help others**

QUART is a proud Supporter of *Medecins sans Frontieres*



[www.quart.de](http://www.quart.de)